

# Pure air generator

## **SPB-3**



### **SPB & GCK series pure air generator character:**

- 1. Two levels pressure stable system, guarantee high precision output pressure
- 2. Magnetic valve automatically automatic changeover, eliminate the noise while water draining out. Water draining and changeover integration.
- 3. Easy operation, water automatically drain out
- 4. Two levels stainless steel made purifications.
- 5. Large inside stainless steel air tank.
- 5. Have two levels pressure stable valve. Output air pressure much stable.
- 6. Have over pressure protection, low voltage startup, DC, greatly protect inside compressor.
- 7. Have A-10 model 0--8L/min & 0--10L/min flow rate pure air generator for choice

## Technical Data:

| Model                     | SPB-3   | GCK3302       | GCK3308       |
|---------------------------|---|---------------|---------------|
| Output Flow               | 0--2L/min   | 0--5L/min     | 0--8L/min     |
| Air Output Purity         | Without oil, two level purification   |               |               |
| Output Pressure (MPa)     | 0--0.4  |               |               |
| Input Power ( W )         | <175  | <300          | <350          |
| Working Condition         | Voltage: 220±20%; 50--60Hz<br>Circumstances Temperature: 15--40 centigrade<br>Relative Humidity: <70%<br>No much dust nor corrosion gas to pollute. |               |               |
| Noisy                     | <35dB   |               |               |
| Pressure Stability        | <0.003Mpa   |               |               |
| Contour Dimension (L*W*H) | 420X250X350mm   | 460X260X350mm | 460*280*350mm |
| Equipment Weight ( kgs )  | <20   | <25           | <25           |