LBK Weighing Scales





Ideal for basic counting and weighing applications, the LBK provides simple operation and a variety of modes. The LBK is easy to use, as it contains just five keys that are colour-coded to highlight frequently used functions. Sleek and compact design makes the LBK ideal for tight work spaces. Solidly constructed with durable ABS plastic housing, the LBK is equipped with a sealed keypad and large stainless steel pan that is easy to remove for cleaning.

The LBK provides the best value for quick, basic weighing tasks







Features

- Vivid, backlit LCD easily visible in any lighting conditions
- Colour-coded keys facilitate quick recognition of the most frequently used buttons
- Level indicator and adjustable feet ensure proper balance setup for optimum weighing results
- · Large, grade 304 stainless steel pan allows easy cleaning
- Lightweight design simplifies transport and portability
- Solid ABS housing is durable for industrial use
- · Audible overload alarm warns when capacity is exceeded
- External calibration allows for verification and adjustment with weights
- Parts counting with preset sample sizes
- · Zero-tracking feature ensures display returns to zero reading
- AC adapter included
- Rechargeable battery included for operation almost anywhere
- · Programmable auto power-off to save energy

Applications

- Weighing
- · Parts counting
- Percentage weighing

About ADAM:

For more than 40 years, professionals worldwide have relied on Adam Equipment for an extensive selection of dependable and affordable weighing equipment.

Adam provides the right balance of speed, performance and value.

LBK Weighing Scales



Model	LBK 3	LBK 6	LBK 12	LBK 30
Capacity	3kg	6kg	12kg	30kg
Readability	0.5g	1g	2g	5g
Repeatability (S.D.)	0.5g	1g	2g	5g
Linearity (+/-)	1g	2g	4g	10g
Pan Size	250x180mm			
Weighing Units	g, kg, lb, lb:oz, oz			
Stabilisation Time (sec)	2			
Calibration	External calibration			
Display	Backlit LCD with 24mm -high digits			
Power Supply	Mains power cord / Internal rechargeable battery			
Operating Temperature	0° to 40°C			
Housing	ABS plastic			
Overall Dim.	250x307x110mm (wxdxh)			
Net Weight	1.8kg			

Accessories

 Item number
 Description

 700200056
 In-use cover

 700200063
 In-use cover (pack of 5)

 700200064
 In-use cover (pack of 10)

 700660290
 Calibration certificate

Key Accessories



In-use cover In-use covers provide an easy and affordable way to help keep the balance clean and protected.

Supplied by