


Hand Held Refractometer MASTER-53α



MASTER-53α are capable of measuring not only various kinds of juice, coffee, and other beverages but also processed food such as sauce, ketchup, low calorie jam, etc. These are an Automatic Temperature Compensation (ATC) model with water resistant feature (IP65). The MASTER-53α is made of metal, and the MASTER-53Pa is plastic. The specifications are exactly the same. For salty and acidic samples, choose the plastic body, MASTER-53Pa.

Specifications

Model	MASTER-53α
Cat.No.	2351
Range	Brix : 0.0 to 53.0%
Minimum Scale	0.2%
Resolution	0.2%
Accuracy	±0.2%
Repeatability	±0.1%

Options

- SUPER EXTRACTOR : RE-29401
- 10% Sucrose (±0.03%) : RE-110010
- 20% Sucrose (±0.03%) : RE-110020
- 30% Sucrose (±0.03%) : RE-110030
- 40% Sucrose (±0.04%) : RE-110040
- 50% Sucrose (±0.05%) : RE-110050
- Daylight plate for small volume samplings : RE-2311-67M
- Daylight plate : RE-2315-60M
- * Calibration Certificate : Contact an ATAGO representative for details.



Shop



Quotation Request



Send an inquiry