

Digital Refractometer PR-101α



This model is suitable for fruits juice, food and beverage as well as chemicals and industrial solutions such as cutting oil, cleaning liquid. The unit is programmed with the user-scale feature which allows to input the coefficient to display the concentrations directly. Also, it is equipped with ELI function for the accurate measurements with confidence even at outside.

Specifications

Model	PR-101α
Cat.No.	3442
Range	Brix : 0.0 to 45.0%
Resolution	0.1%
Accuracy	±0.1%
Temperature compensation	5 to 40°C
Power Supply	006P dry battery (9V)
Dimensions & Weight	17×9×4cm, 270g



Shop



Quotation Request



Send an inquiry

Options

- SUPER EXTRACTOR : RE-29401
- 10% Sucrose (±0.03%) : RE-110010
- 20% Sucrose (±0.03%) : RE-110020
- 30% Sucrose (±0.03%) : RE-110030
- 40% Sucrose (±0.04%) : RE-110040

* Calibration Certificate : Contact an ATAGO representative for details.