

Digital Refractometer PR-201α



It is suitable for various samples up to Brix 60% including jams, concentrated juices and lubricant. The unit is programmed with the user-scale feature which allows to input the coefficient to display the concentrations directly. Also, it is equipped with ELI function for the accurate measurements with confidence even at outside.

Specifications

Model	PR-201α
Cat.No.	3452
Range	Brix : 0.0 to 60.0%
Resolution	0.1%
Accuracy	±0.1%
Temperature compensation	10 to 40°C
Power Supply	006P dry battery (9V)
Dimensions & Weight	17×9×4cm, 270g



Shop



Quotation Request



Send an inquiry

Options

- 10% Sucrose (±0.03%) : RE-110010
- 20% Sucrose (±0.03%) : RE-110020
- 30% Sucrose (±0.03%) : RE-110030
- 40% Sucrose (±0.04%) : RE-110040
- 50% Sucrose (±0.05%) : RE-110050
- 60% Sucrose (±0.05%) : RE-110060

* Calibration Certificate : Contact an ATAGO representative for details.