

## Digital Refractometer PR-301α



Shop



Quotation Request



Send an inquiry

It is suitable for higher concentration like jellies, raw liquid sugar, etc. The unit is programmed with the user-scale feature which allows to input the coefficient to display the concentrations directly. Also, it is equipped with ELI function for the accurate measurements with confidence even at outside.

#### Specifications

Model	PR-301α
Cat.No.	3462
Range	Brix : 45.0 to 90.0%
Resolution	0.1%
Accuracy	±0.1%
Temperature compensation	10 to 40°C
Power Supply	006P dry battery (9V)
Dimensions & Weight	17×9×4cm, 270g
Included Accessories	Standard liquid LG

#### Options

- 50% Sucrose (±0.05%) : RE-110050
- 60% Sucrose (±0.05%) : RE-110060
- \* Calibration Certificate : Contact an ATAGO representative for details.

#### Parts

- Standard liquid LG (7mL) with round sheet : RE-9325