

## Digital Butyro Refractometer PR-BUTYRO



New PR-BUTYRO perfect for Butter, Ghee & Edible Oil industry for checking purity / adulteration. PR-BUTYRO has dual scales of Butyro and Refractive Index. The temperature is automatically corrected and not necessary to worry about the temperature.

#### Specifications

Model	PR-BUTYRO
Cat.No.	3454
Range	Butyro : 30.0 to 90.0 RI : 1.4450 to 1.4850 (Converted at 40°C)
Resolution	Butyro: 0.1 RI : 0.0001
Accuracy	Butyro : $\pm 0.5$ (at 40°C) RI : $\pm 0.0003$ (at 40°C)
Temperature compensation	10 to 50°C
Power Supply	006P dry battery (9V)
Dimensions & Weight	17×9×4cm, 270g


[Shop](#)

[Quotation Request](#)

[Send an inquiry](#)

#### Options

\* Calibration Certificate : Contact an ATAGO representative for details.