

Digital Refractometer PR-RI



It is suitable for measuring the refractive indices of liquid samples in pharmaceutical and chemical products. When the start key is pressed after placing some droplets of sample on the prism, the refractive index and temperature will be displayed.

Specifications

| | |
|---------------------|--|
| Model | PR-RI |
| Cat.No. | 3480 |
| Range | Refractive index : 1.3306 to 1.4436 Temperature 5 to 45°C |
| Resolution | 0.0001 |
| Accuracy | ±0.0002 (water at 20°C) |
| Ambient Temperature | 10 to 40°C |
| Power Supply | 006P dry battery (9V) |
| Dimensions & Weight | 17×9×4cm, 270g |



Shop



Quotation Request



Send an inquiry

Options

* Calibration Certificate : Contact an ATAGO representative for details.