

Digital Urine S. G. Refractometer UG- $\alpha$ 

This new Digital Urine Specific Gravity Refractometer UG- $\alpha$  (alpha) is an upgraded version of the late UG-1. After pressing the START key, with the sample placed on the prism, the measurement values will be displayed within 3 seconds.

#### Specifications

Model	UG- $\alpha$
Cat.No.	3464
Range	Urine S.G. : 1.0000 to 1.0600
Resolution	0.0001
Accuracy	$\pm 0.0010$
Temperature compensation	10 to 35°C
Power Supply	006P Dry battery (9V)
Dimensions & Weight	17×9×4cm, 270g



Shop



Quotation Request



Send an inquiry

#### Options

\* Calibration Certificate : Contact an ATAGO representative for details.